

KOMENDA	OPIS ŻĄDANIA	UWAGI
0x4D	CMD_SUSPEND_RESUME	ACK
0x5D	CMD_POWER_CTRL	ACK
0x70	CMD_READ_RTC	
0x71	CMD_READ_PUMP_ID	
0x72	CMD_READ_BATTERY	
0x73	CMD_READ_RESERVOIR	
0x74	CMD_READ_FIRMWARE_VERSION	
0x75	CMD_READ_ERROR_STATUS	
0x76	CMD_READ_REMOTE_CTRL_IDS	
0x79	CMD_READ_TOTALS_TODAY	
0x80	CMD_READ_HISTORY_DATA	ACK
0x81	???	
0x83	CMD_READ_PUMP_STATE	
0x86	CMD_READ_LANGUAGE	
0x87	CMD_READ_BOLUS_WIZARD_SETUP_STATUS	
0x88	CMD_READ_CARB_UNITS	
0x8A	CMD_READ_CARB_RATIOS	
0x8B	CMD_READ_INSULIN_SENSITIVITIES	
0x8D	CMD_READ_PUMP_MODEL_NUMBER	
0x8E	CMD_READ_BG_ALARM_CLOCKS	
0x8F	CMD_READ_RESERVOIR_WARNING	
0x90	CMD_READ_BG_REMINDER_ENABLE	
0x92	CMD_READ_STD_PROFILES	
0x93	CMD_READ_A_PROFILES	

KOMENDA	OPIS ŻĄDANIA	UWAGI
0x94	CMD_READ_B_PROFILES	
0x95	CMD_READ_LOGIC_LINK_IDS	
0x97	CMD_READ_BG_ALARM_ENABLE	
0x98	CMD_READ_TEMP_BASAL	
0x9A	CMD_READ_GLUCOSE_HISTORY	ACK
0x9B	CMD_READ_ISIG_HISTORY	ACK
0x9C	CMD_READ_CALIBRATION_FACTOR	
0x9D	CMD_READ_CURRENT_HISTORY_PAGE_NUMBER	
0x9E	???	
0x9F	CMD_READ_BG_TARGETS	
0xA0	???	ACK
0xA3	CMD_READ_PUMP_TRACE	
0xA4	CMD_READ_DETAIL_TRACE	
0xA6	CMD_READ_NEW_TRACE_ALARM	
0xA7	CMD_READ_OLD_TRACE_ALARM	
0xA8	CMD_FILTER_BG	ACK
0xA9	CMD_FILTER_ISIG	ACK
0xAB	???	
0xAC	???	
0xAD	???	ACK
0xC0	CMD_READ_SETTINGS	
0xC1	CMD_READ_SAVED_SETTINGS_DATE	
0xC2	???	
0xC3	CMD_READ_CONTRAST	

KOMENDA	OPIS ŻĄDANIA	UWAGI
0xC5	CMD_READ_BOLUS_REMINDER_ENABLE	
0xC6	CMD_READ_BOLUS_REMINDERS	
0xC7	CMD_READ_FACTORY_PARAMETERS	
0xC8	???	
0xC9	???	
0xCD	CMD_READ_CURRENT_GLUCOSE_HISTORY_PAGE_NUMBER	
0xCE	CMD_READ_CURRENT_PUMP_STATUS	
0xCF	CMD_READ_SENSOR_SETTINGS	
0xD1	CMD_READ_SENSOR_PREDICTIVE_ALERTS	
0xD2	CMD_READ_SENSOR_DEMO_AND_GRAPH_TIMEOUT	
0xD3	CMD_READ_SENSOR_ALARM_SILENCE	
0xD4	CMD_READ_SENSOR_RATE_OF_CHANGE_ALERTS	
0xD5	CMD_READ_VCNTR_HISTORY	ACK
0xDB	???	
0xDD	???	
0xDF	???	
0xE6	???	
0xE8	???	ACK
0xEC	???	
0xEF	???	
0xF0	CMD_READ_OTHER_DEVICES_IDS	
0xF1	???	
0xF3	???	